mesofiller® periocular

INSTRUCTIONS

DESCRIPTION

mesofiller® is a sterile, injectable, non-pyrogenic, reabsorbable medical product made out of cross-linked hyaluronic acid of non-animal origin, produced via bacterial fermentation. **mesofiller®** is a colourless gel contained in a pre-filled, graduated, disposable, and sterile syringe with a Luer Lock adapter.

COMPOSITION

mesofiller® periocular

- Cross-linked hyaluronic acid 15 mg/ml
- Saline sterile solution as needed

MECHANISM OF ACTION

mesofiller®'s action consists in increasing the volume of the dermo-epidermal tissue based on the natural ability of the hydrophile molecules of hyaluronic acid to bind to quantities of water that are many times heavier than they are. This allows them to fill the intradermal spaces and fill the intercellular matrix, making the tissues firmer. The reticulation of the hyaluronic acid contained in the product makes the filler effect more stable and lasting over time.

RECOMMENDATIONS

mesofiller® is intended to be used as a temporary filler to correct small skin defects, such as wrinkles or scars. It is also intended for the treatment of pathological conditions such as HIV-associated lipodystrophy.

In particular, the following use are recommended:

- **mesofiller® periocular** for the treatment of periocular wrinkles.

CONTRAINDICATIONS

mesofiller® must not be used:

- for injections in the palpebral periorbital area;
- for injections into blood vessels;
- for breast augmentation, or for bone, tendon, ligament or muscular implantation.

mesofiller® must not be administered to patients with:

- ascertained tendency to develop hypertrophic or keloid scars or pigmentation disorders;
- previous active autoimmune diseases or those under immunotherapy;
- serious multiple allergies or active inflammations;
- coagulation disorders;
- acknowledged hypersensitivity to hyaluronic acid.

mesofiller® must not be used on pregnant or breast feeding women nor on people below the age of eighteen.

MEDICAL WARNING

mesofiller® must only be used for subcutaneous or intradermal injection and must not be injected intravasally. There is the potential risk of the material being accidentally injected into blood vessels when executing the procedure. In rare cases, this may cause vascular occlusion, abscess, necrosis or embolic phenomena. In order to verify that the needle is not in a blood vessel during the normal use of the device, lightly aspirate with the syringe before proceeding with the injection.

mesofiller® is packaged to be used on only one patient and within the same session.

mesofiller® is provided as a sterile, disposable product and must not be resterilised. Do not reuse.

Use **mesofiller®** straight after opening the blister and dispose of the syringe immediately after use, even if the entire content has not been used.

Do not use the product after the expiry date printed on the package.

Do not use **mesofiller®** with other injectable products and do not mix it with other preparations before use.

Do not use **mesofiller®** when the skin presents reddened areas, swelling or pigment alteration; do not use **mesofiller®** during

inflammatory, infectious or feverish states.

Patients that use drugs to slow down blood coagulation (such as aspirin) may present ekhymosis or pronounced bleeding around the site of the injection. Patients being treated with anticoagulants shall have to refrain from taking them, in agreement with their physician, at least 14 days before receiving **mesofiller®** injections.

Keep the product out of reach from children.

PRECAUTIONS

Use **mesofiller®** only according to these instructions.

mesofiller® is a colourless transparent gel: if the content of the syringe appears cloudy or contains particles, do not use the product and promptly notify **mesoestetic Pharma Group, S.L** at the address provided below.

Do not use the device if the sterile syringe or the blister containing it appear to be damaged.

Do not use the needles provided if they are bent and/or deformed or if the packaging is not intact.

Do not use non-sterile needles.

Take all normal precautions associated with intradermal injection procedures. The injection is intrinsically associated to risk of infection. Skin must be perfectly cleansed and disinfected before treatment. Do not use on patients with active infections or inflammatory processes in the vicinity or on the injection delivery site.

Do not use **mesofiller®** in areas where other implantations have been carried out and do not use at the same time as any other injectable product.

Advise patients not to wear makeup during the 24 hours following the injection and inform them that, during the remission period, meaning until swelling and redness has completely disappeared, they must avoid exposing the treated area to direct sunlight, UV lamps or to the cold.

Inflammatory reactions may occur if laser, chemical peels, or dermo-abrasive treatments are given after using **mesofiller®**. This is why **mesofiller®** should only be used when the effects of any prior treatment have disappeared completely.

Syringes and needles may be a potential biological hazard once they have been used. Handle them with care and dispose of them in compliance with approved medical practices.

NEGATIVE REACTIONS

Very important: before beginning treatment, the doctor must inform the patient that there are potential side effects linked to intradermal implantation that may occur immediately or some time following treatment.

Following the use of **mesofiller®**, some side effects may occur related to the injection, including pain, itching, alterations in pigmentation and density of the tissue in the relevant injection site, which normally disappear spontaneously within one week. Small oedemas may sometimes appear, but should reabsorb within a few days. Among the negative reactions that are correlated to the use of the product, according to bibliographical studies of the existing literature and the scientific knowledge of the manufacturer as at the date of this information leaflet, we point out the following:

- inflammatory reactions, itching, painful to the touch following injection, lasting up to a maximum of one week;
- bleeding or haematoma in rare cases;
- follicular papules-pustules that may form up to four weeks following treatment lasting for two weeks;
- allergic reactions to hyaluronic acid;
- hardening or appearance of nodules at the injection site;
- in extremely rare cases, granulomas or significant inflammatory reactions, wheals, necrosis, rash, which may even appear with a 2-4 week delay.

IMPORTANT: patients must report any type of negative reaction, whether serious or not, listed in this information leaflet, and the doctor, as well as prescribing appropriate treatment to cure such reactions, is under the obligation to promptly notify the distributor and manufacturer, (**mesoestetic Pharma Group, S.L**).

INCOMPATIBILITY: there is an acknowledged incompatibility between hyaluronic acid and quaternary ammonium salts, such

as benzalkonium chloride. Therefore, avoid placing **mesofiller®** in contact with such substances, or using medical and surgical equipment containing such quaternary ammonium salts.

PROCEDURE

Very important: before beginning treatment, the doctor must inform the patient on the composition, characteristics, precautions, medical warnings, contraindications and potential negative effects associated to this device. Furthermore, the doctor must take a medical history and assess any prior immunological or allergic problems of the patient before proceeding, to ascertain compatibility with the use of **mesofiller®**. Finally, the doctor must assess the relevant areas of the skin to be treated to ensure that the use of such treatment is appropriate and useful for the patient.

mesofiller® may only be administered by qualified medical personnel. For a correct final result of treatment, it is essential to proceed with a correct injection technique.

Before the injection, disinfect the area with alcohol or other antiseptic solutions.

The product must be injected using a sterile syringe with a thin needle. We recommend using a needle cannula 25G 2" or a needle 30G 1/2" for **mesofiller® periorcular**. The choice and use of the needle and injection procedure is in any case left to the doctor's discretion.

The quantity of the product to be injected must be established according to the depth of the wrinkle, and may be easily dosed thanks to the use of graduated needles.

Straight after injection, a firm pressure and massage must be applied to the treated area so as to favour a smooth distribution of the product.

Before each procedure, one must attentively assess the symmetry of the facial area, planning to inject the same quantity of material in the right part as in the left, and in symmetrical points. Immediately before proceeding with the implantation, remove the syringe from its protective blister, unscrew the cap of the syringe and place the needle onto the Luer Lock adapter by screwing on the base clockwise until it locks.

To ensure that **mesofiller®** is used safely, it is important to use an appropriate needle and to assemble the needle and syringe correctly, as an incorrect assembly may cause the needle to separate from the syringe during the injection (see figures 1 and 2).


Fig. 1
Unscrew the cap of
the Luer Lock adapter


Fig. 2
Insert the needle

Remove the protective cap of the needle only just before making the injection.

Before injecting, expel the air from the syringe until a drop of the product appears on the tip of the needle.

Introduce the needle with a 30 degree angle, following the exact direction of the wrinkle or groove, and inject the product into the mid derma. If the hyaluronic acid is injected in too much depth, the final effect will last less and, on the contrary, if it is injected too close to the surface, nodules and granulomas may form.

Never apply too much pressure on the syringe; inject **mesofiller®** whilst slowly retracting the needle, depositing a homogenous quantity of the product along the chosen site. Interrupt the pressure on the plunger of the syringe before extracting the needle from the skin, in order to prevent the product from spilling out of the point of injection. In case of immediate swelling it is useful to briefly apply cold packs to the area.

DOSAGE

Dosage depends on the skin condition of the patient and depth of the treated area. Normally 0.5-1.0 ml per site is necessary. Do not exceed

the quantity of one syringe (1 ml) per site and per session. Optimal results may be obtained by repeating the treatment regularly over time, considering the various factors that contribute to the duration of the effects of the injections, such as life-style, physical activity, sunlight exposure or tanning lamps, age and appropriate filling technique.

PACKAGING

One package of **mesofiller® periocular** contains:

- One pre-filled, graduated, 1 ml syringe containing the product with a Luer Lock adapter sealed in a blister.
- 1 sterile disposable needle cannula 25G 2" (CE 0123).
- 1 sterile disposable needle 30G 1/2" (CE 0123).
- Plastic support assembled with the syringe in order to make it easier to handle.
- Set of adhesive labels with the lot number, one of which must be applied to the patient's clinical file, and the other delivered to the patient. These procedures are strictly compulsory in order to ensure the traceability of the product.
- One information leaflet that complies with European regulations to be read attentively before using the device.

STERILISATION

The content of **mesofiller®** syringes is steam sterilised.

The needles are sterilised by irradiation.

STORAGE CONDITIONS

Read the expiry date of the product on the syringe, on the label of the blister or on the package. Store in a dry environment, shielded from light and at a temperature between 5°C and 25°C. Fragile.

Note: the graduation on the syringe is not for measurement purposes but is only provided as a visual control over the injected material.

SYMBOLS ON PACKAGE

	Do not reuse, one use only
	Follow instructions attentively
	Expiry date
	Lot number
	Steam sterilization
	Sterilized using irradiation
	Store in dry environment
	Shield from sunlight
	Store at a temperature between 5° and 25° C
	Manufacturer
	Date of manufacture
	Product code
	Do not use if packaging is damaged
	Fragile, handle with care
	Syringe
	Needle
	CE marking in compliance with the 93/42/EEC directive concerning medical devices and subsequent amendments and integrations. 0373 is the number of the Notified Body.

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